
# FlexiFit<sup>14</sup>31

FULL FACE MASK



### **Specification Sheet**

A Full Face Delivery Mask






Incorporating the FlexiFit™
Technology with auto-contouring
properties, the FlexiFit™ 431
Full Face Mask comes with
an Under-Chin fit and highly
contoured Silicone Seal for
maximum comfort and stability.





## Full Face Delivery - *FlexiFit*™431




The FlexFit<sup>™</sup> 431 Full Face Mask offers a comfortable alternative to nasal or oral delivery and comes as an Out-of-Box Solution fully assembled and ready for use.

#### UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



#### FlexiFit<sup>™</sup> Technology Auto-contouring fit

With Air Cushion
FlexiFit™ Technology
automatically adjusts to a
wide range of nasal and facial
contours with simple onestep ease of fitting and no
complicated adjustments.



#### Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.




## Jnder-Chin Design

An Under-Chin fit provides stability by providing a constant point of reference and maintains the seal even with movement of the jaw during the night.



## Efficient Bias Diffuser Ouiet usage

Maintenance Free
This advanced air diffuser
system prevents air from
blowing onto the sleeping
partner and decreases noise that
could potentially disrupt sleep.



## Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.



## Crown Strap Extra stability

Complements the Under-Chin Design to maximize stability and eliminate the need for a chin strap.



#### Easy-Release Strap One-step removal

Simply pull the Easy-Release Cord for one-step removal of headgear. No ongoing headgear adjustments required.

#### CONTENTS

- Complete mask with headgear assembled and ready for use
- Alternate size Silicone Seals included (S,M,L)
- User friendly Fitting Guide

REORDER CODE: HC431

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER AND PAYKEL HEALTHCARE REPRESENTATIVE: