



For more information, please contact your local Fisher & Paykel Healthcare representative

Performs Automatically - in Tune with You

Manufacturer
Fisher & Paykel Healthcare Ltd
15 Maurice Paykel Place
East Tamaki, Auckland 2013
PO Box 14 348 Panmure
Auckland 1741
New Zealand
Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Australia
Fisher & Paykel Healthcare Pty Limited
36-40 New Street, PO Box 167
Ringwood, Melbourne
Victoria 3134, Australia
Tel: +61 3 9879 5022
Fax: +61 3 9879 5232

Austria
Tel: 0800 29 31 23
Fax: 0800 29 31 22

Benelux
Tel: +31 40 216 3555
Fax: +31 40 216 3554

China
Tel: +86 20 3205 3486
Fax: +86 20 3205 2132

France
Tel: +33 1 6446 5201
Fax: +33 1 6446 5221

Germany
Tel: +49 7181 98599 0
Fax: +49 7181 98599 66

India
Tel: +91 80 4284 4000
Fax: +91 80 4123 6044

Irish Republic
Tel: 1800 409 011

Italy
Tel: +39 06 7839 2939
Fax: +39 06 7814 7709

Japan
Tel: +81 3 3661 7205
Fax: +81 3 3661 7206

Northern Ireland
Tel: 0800 132 189

Spain
Tel: +34 902 013 346
Fax: +34 902 013 379

Sweden
Tel: +46 8 564 76 680
Fax: +46 8 36 63 10

Switzerland
Tel: 0800 83 47 63
Fax: 0800 83 47 54


Taiwan
Tel: +886 2 8751 1739
Fax: +886 2 8751 5625

Turkey
Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK **EC** **REP**
Fisher & Paykel Healthcare Ltd
Unit 16, Cordwallis Park
Clivemont Road,
Maidenhead
Berkshire SL6 7BU, UK
Tel: +44 1628 626 136
Fax: +44 1628 626 146

USA/CANADA
Tel: +1 800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001

フィッシャー & パイケルヘルスケア株式会社
〒103-0025 東京都中央区日本橋茅場町3-8-10 リベラ茅場町5F
FAX: 03-3661-7206 注文専用 FAX: 0800-100-7207
TEL: 03-3661-7205 (代表)



F&P ESON NASAL MASK

F&P Eson is designed to automatically perform in tune with you and your needs.

Three small, streamlined mask components work in harmony for greater performance, stability and auto-adjusting comfort. This compact nasal mask is packed with technological innovations orchestrated to make SleepLife easier.



Mask Reorder Codes		Spare Parts Reorder Codes			
400449	F&P Eson Mask Small	400HC567	F&P Eson Headgear - Small	400HC228	F&P Eson Diffuser (10pk) and Cover
400450	F&P Eson Mask Medium	400HC568	F&P Eson Headgear - Medium/Large	400HC570	F&P Eson Swivel
400451	F&P Eson Mask Large	400HC569	F&P Eson Headgear - Clips & Buckle	400HC571	F&P Eson Mask No H/G - Small
		400HC126	F&P Eson Seal - Small	400HC572	F&P Eson Mask No H/G - Medium
		400HC127	F&P Eson Seal - Medium	400HC573	F&P Eson Mask No H/G - Large
		400HC128	F&P Eson Seal - Large		

F&P Eson is very 'easy on you' from bag to bed. It has been designed to perform automatically in tune with every aspect of SleepLife.* F&P Eson maximizes no-fuss comfort and simplicity by minimizing manual adjustments, parts and noise.

EASY ON PROVIDERS AND SLEEPLABS

EASY ON PATIENTS AND THEIR PARTNERS

Easy Fit



Intuitive fitting innovations make it easy and fast to fit your patients.



Easy-Clip hooks allow for fast fitting and removal.

Easy Clean



With only three simple parts, it is easy to clean and reassemble.

Easy Comfort



Design innovations maximize comfort and optimize a stable, effective seal all night.



Easy Sleep



It is super-quiet with minimal draft allowing an undisturbed sleep for you and your partner.



F&P Eson's three simple components, the **RollFit™ Seal**, **ErgoFit Headgear** and **Easy Frame**, work in harmony to deliver the comfort, seal and easy use that Fisher & Paykel Healthcare masks are known for.



1. RollFit Seal

Easy auto-adjusting fit that optimizes an effective, comfortable seal

As the name suggests, the one-piece seal 'rolls' back and forth on the bridge of the nose to adjust automatically. In doing so, the **RollFit Technology** minimizes pressure on the bridge of the nose.

The **RollFit Seal** incorporates **Integrated Stability Panels** for optimal sealing performance and stability all night long.

2. ErgoFit Headgear

Easy to use, fit, adjust and clean

This breathable, ergonomically designed headgear has been designed to self-locate high on the rear of the head allowing for maximum head movement (sideways, up and down) without mask dislodgement.

The **ErgoFit Headgear** also incorporates **Embedded VELCRO®* Tabs (2a)** with **Easy-Clip hooks (2b)** for fast and easy fitting and removal. The adjustable **Crown Strap (2c)** enhances mask stability.

* VELCRO® is a registered trademark of Velcro Industries B.V.

3. Easy Frame

Easy to wear, see and sleep

This low-profile frame is stable, durable and small, and ensures a clear line of sight.

The one frame fits all three seal sizes and has an **'Easy-Clip'** frame attachment for effortless assembly after cleaning.

The **Easy Frame** comes with the **Q Cover (3a)** for super-quiet use and minimal partner draft.

The **Ball and Socket Elbow** rotates for freedom of movement while minimizing drag from the CPAP tube.

F&P Eson performs in tune with what patients and their healthcare providers want: **Automatic Performance** that makes adherence **Easy on Everyone.**