



DeVilbiss IntelliPAP® Standard CPAP System

You talked. We listened.

Everyone gets what they need.

The DeVilbiss IntelliPAP Standard CPAP System was designed with patients and providers in mind. IntelliPAP's combination of comfort, education and usage tracking with SmartCode® helps to ensure patient compliance.

From its reliable electronics and water-resistant housing to its durable Lexan® humidifier chamber and stainless steel heat transfer plate, IntelliPAP is built for years of dependable performance and patient satisfaction. It's just what you would expect from DeVilbiss.



DV51D Series

Performance

- **Quality** – Less than 1% failure rate and backed by the best-in-industry, 3-year warranty
- **Compliance** – Free remote tracking of compliance data through SmartCode technology
- **Comfort** – Optimal heated humidification system
- **Design** – Stacked design creates an exceptionally small footprint and whisper-quiet operation [26 dBA]*

DeVilbiss IntelliPAP Standard CPAP System

Design

- Easy and intuitive setup
- Auto on/off** automatically starts flow of air with patient breath
- Eight-foot power cord provides greater positioning flexibility – no bulky power brick
- Reminders alert patient to replace or check parts
- Includes attractive travel bag
- Large, slip-resistant pads on bottom of unit to minimize movement on nightstand



Humidifier Chamber

- Overfill protection
- Corrosive resistant, stainless steel plate
- Unique baffling increases air water exchange for optimal humidity
- 400 ml capacity for 8 hours of use at the highest setting



Usage Data	
	Last 90 Days
% Days of at least 4 hours (not continuous) in a 24 hour period	54%
Day Count	30
Days at least 4 hours	58
SmartCode	2MTJT-45G4-7DVG8
White Breathing Hours Last Day	
Adherence Score (4 continuous hours)	67% (28/33)

SmartCode

- Encrypted usage data codes that can be communicated over the phone, e-mail, or text to generate compliance and clinical reports



NEW EasyFit® SilkGel Mask is the perfect fit with IntelliPAP



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com

* Tested per ISO 17510 ** U.S. Patent 5,503,146