

RESMED




RESMED

Mirage Micro™

NASAL MASK

Micro makes (f)it easy

Visit www.myresmed.com for more information about the Mirage Micro and other ResMed products, as well as fitting guides and other therapy tips.



ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Belgium, Brazil, China, Finland, France, Germany, Greece, Hong Kong, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Norway, Portugal, Singapore, South Africa, Spain, Sweden, Switzerland, United Kingdom (see website for details).

Mirage Micro is a trademark of ResMed Ltd. and Mirage is registered in the U.S. Patent and Trademark Office. ©2009 ResMed. 1012294/2 09 01

The Mirage Micro™ nasal mask combines personalized fit and easy-to-use features to ensure your therapy experience is comfortable and effective.

Personalized fit and comfort

- **Micro makes (f)it easy** – MicroFit dial allows for up to 24 precise adjustments of the forehead support, setting the mask at the most comfortable position.
- **Proven seal** – Mirage™ dual-wall cushion technology reduces pressure on the nasal bridge and seals effectively.
- **A winning combination** – As you adjust the dial to your profile, the cushion works to seal effectively, so there is no need to over-tighten the headgear.
- **A range of sizes** – With five cushions to choose from, you'll find the perfect fit for you.

Quiet, sleek and unobtrusive

- **Quiet performance** – Unique radial design quietly and gently disperses air 360°.
- **Clear field of vision** – Slender forehead support is unobtrusive and feels light on the face.
- **Designed for your lifestyle** – Softly diffuses airflow so you and your partner can enjoy a good night's sleep.
- **Includes convenient drawstring bag** to fit your lifestyle and travel needs.



Intuitively easy

- **First-time fit** – It's easy to adjust to therapy when your mask fits you well the first time you put it on. Your Mirage Micro is designed to fit comfortably the first time and every time.
- **Ready to wear** – Handy headgear clips are easy to attach and detach, saving your optimal headgear settings so you can "set and forget!"
- **One-handed adjustment** – There's no need to take off your mask to use the MicroFit dial. You can simply adjust the dial with one hand while wearing the mask, and even while lying down.
- **Easy to care for your mask** – Fewer and simplified parts make your mask easy to assemble, disassemble and clean.



Amazingly quiet and diffused new vent design