

Swift™ FX Nano... small in size,
big on comfort.



Swift FX Nano is the new easy breathing compact mask your patients have been waiting for! It comes with a soft, low-profile nasal cushion and minimal headgear that really opens up the face. It's time for your patients to slip into something more comfortable.

Benefits for you ...

Builds on the strength of Swift FX

Swift FX Nano delivers the compact size, facial freedom, reliability and comfort that the Swift name is famous for with the breathing ease of a nasal cushion.

Set up once and they're set

It fits quickly and requires minimal or often no adjustment once fitted, saving time, money and additional patient visits.

Simple to use and reassemble after cleaning

Swift FX Nano has fewer main parts than most other nasal masks (just three: headgear, cushion and tube) and it is simple to use.

Maintains a comfortable seal

Its smooth 360° rotating ball-joint helps maintain seal effectively and provides comfort with adaptive flexibility.

Designed for easy user acceptance

Swift FX Nano is light and compact with minimal headgear helping encourage patients to accept and continue therapy.

... and your patients

Light and natural

It provides the lightness of a Swift FX pillows mask with the natural breathing comfort of a nasal mask.

Can easily wear glasses on therapy

Swift FX Nano provides uninterrupted vision, so patients can wear glasses to comfortably watch TV, or read, while on therapy.

Choice of sleeping positions

Swift FX Nano is flexible during movement as the user sleeps—so patients can sleep in a variety of different positions.

Surprisingly quiet performance

Advanced venting design diffuses airflow so the mask is quiet to use—just 25 dBA.*

* Declared dual numbered noise emission values (with uncertainty 3 dBA) in accordance with ISO 4871.e1.



Swift™ FX Nano
NASAL MASK

Swift™ FX Nano for Her
NASAL MASK



Small, light and appealing

With its compact, open-faced design, minimal 2-point headgear (no need for rigid frames or forehead support) light weight, and optional soft wraps, Swift FX Nano provides patients with comfort and a real sense of facial freedom.

Sleep apnea patients often look for the smallest, least intrusive and most comfortable mask they can find. With the launch of Swift FX Nano, patients can now enjoy the natural breathing comfort of a nasal mask in an appealing new compact size.

Quick fitting headgear
with easy adjusting notches on the top strap. Set once and forget!

SoftEdge™ fabric back strap
is soft and comfortable, with quick setting fasteners

All soft nasal cushion
with a dual-wall design and Spring-Air Technology™ that is soft for the patient while providing a comfortable seal

Ultra smooth ball-joint
with 360° unimpeded movement for adaptive flexibility

Whisper-quiet diffuse venting
gently distribute air away from the user and their bed partner

Innovative cuff
decouples the turning force from elbow to the short tube, allowing the tube to move independently of the elbow

Product Codes

	Europe 1	Europe 2	Europe 3
Mask system* (Standard)	62202	62204	62206
Mask system* (Wide)	62253	62255	62257
For Her mask system** (Small)	62203	62205	62207

* Includes Standard or Wide cushion, Swift FX Nano headgear assembly, plus grey backstrap, soft wraps and frame assembly.
** Includes Small cushion, Swift FX Nano For Her headgear assembly, plus pink backstrap, soft wraps and frame assembly.

Europe 1: English (UK), French, German and Italian (Switzerland)

Europe 2: English, Danish, Finnish, Norwegian and Swedish

Europe 3: English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish

ResMed Ltd
1 Elizabeth Macarthur Drive
Bella Vista NSW 2153 Australia

Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA. **EC/REP** ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See www.resmed.com for other ResMed locations worldwide. Swift, SoftEdge and Spring Air Technology are trademarks of ResMed Ltd and are registered in the U.S. Patent and Trademark Office. © 2013 ResMed Ltd. Specifications may change without notice. 1015961/1 2013-01



2-point headgear
has only minimal contact on the face and provides excellent all round vision



'Lock and key'
for easy attaching and detaching of headgear

Size and name indicators
clearly identify the size and the correct orientation of the mask

Cushion tab hinges
adapt to improve flexibility between headgear and cushion

Spring-flex tubing
minimises any pull on the mask, enhancing comfort on the face

Swift™ FX Nano
NASAL MASK

Swift™ FX Nano for Her
NASAL MASK