





Ultra Mirage™ II
NASAL MASK


Mirage Activa™
NASAL MASK




↑
STANDARD



↑
LARGE



↑
SHALLOW



↑
SHALLOW WIDE
(Ultra Mirage II only)



Actual Size: 5 1/2" (Width) x 7" (Height)
Ruler is for scaling purposes.

- How to use the fitting template:
1. Cut out the template along the dotted line
 2. Cut out nasal holes
(hint: fold each nasal hole in half to cut)
 3. Follow directions on template to determine appropriate cushion size

RESMED